


Product Datasheet

Superoxide Dismutase 3/EC-SOD/Sod3 Antibody (orb763080)


Description

Superoxide Dismutase 3/EC-SOD/Sod3 Antibody


Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Immunogen	E.coli-derived mouse Sod3 recombinant protein (Position: V36-L101).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	32 kDa
Uniprot ID	O09164
Dilution Range	Western blot, 0.25-0.5µg/ml, Mouse Immunohistochemistry (Paraffin-embedded Section), 2-5µg/ml, Mouse Flow Cytometry, 1-3µg/1x10 ⁶ cells, Mouse Direct ELISA, 0.1-0.5µg/ml, Mouse
Expiration Date	12 months from date of receipt.



Western blot analysis of Superoxide Dism...



IHC analysis of Superoxide Dismutase 3/E...



Flow Cytometry analysis of HEPA1-6 cells...