



Product Datasheet

RPL13A Antibody (monoclonal, 4C9) (orb763102)


Description	RPL13A Antibody (monoclonal, 4C9)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, WB
Immunogen	E.coli-derived human RPL13A recombinant protein (Position:A2-V203).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	24 kDa
Uniprot ID	P40429
Dilution Range	Western blot, 0.25-0.5µg/ml, Human Immunocytochemistry/Immunofluorescence, 5µg/ml, Human Flow Cytometry, 1-3µg/1x10 ⁶ cells, Human
Expiration Date	12 months from date of receipt.



Western blot analysis of RPL13A using an...



IF analysis of RPL13A using anti-RPL13A



Flow Cytometry analysis of A431 cells us...