


Product Datasheet


IL-2 Receptor alpha/Il2ra Antibody (orb763177)

Description
IL-2 Receptor alpha/IL2ra Antibody


Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence of mouse IL-2 Receptor alpha/IL2ra (PWKHEDSKRIYHFVE).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	50 kDa
Uniprot ID	P01590
Dilution Range	Western blot, 0.25-0.5 µg/ml, Mouse, Rat Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Mouse, Rat Flow Cytometry, 1-3 µg/1x10 ⁶ cells, Mouse
Expiration Date	12 months from date of receipt.



Western blot analysis of IL-2 Receptor a...



IHC analysis of IL-2 Receptor alpha/IL2r...



IHC analysis of IL-2 Receptor alpha/IL2r...