

Product Datasheet

TBP-1/PSMC3 Antibody (monoclonal, 4D3) (orb763201)

Description

TBP-1/PSMC3 Antibody (monoclonal, 4D3)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, ICC, IF, IHC, WB

Immunogen

E.coli-derived human TBP-1/PSMC3 recombinant protein (Position: M1-A439).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW

50 kDa

Uniprot ID


P17980

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat
Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human
Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



Western blot analysis of TBP-1/PSMC3 using anti-TBP-1/PSMC3 antibody.



IHC analysis of TBP-1/PSMC3 using anti-TBP-1/PSMC3 antibody.



IHC analysis of TBP-1/PSMC3 using anti-TBP-1/PSMC3 antibody.