

Product Datasheet

RPA32/RPA2 Antibody (monoclonal, 3B2E9) (orb763203)

Description

RPA32/RPA2 Antibody (monoclonal, 3B2E9)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, ICC, IF, IHC, WB

Applications
Immunogen

E.coli-derived human RPA32/RPA2 recombinant protein
(Position: Q34-H254).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of
500 µg/ml.

Storage

At -20°C for one year from date of receipt. After
reconstitution, at 4°C for one month. It can also be
aliquotted and stored frozen at -20°C for six months. Avoid
repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive
results. Other applications have not been tested. Optimal
dilutions should be determined by end users. Adding 0.2
ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2a

Clonality

Monoclonal

MW


32 kDa

Uniprot ID
P15927
Dilution Range


Western blot, 0.25-0.5 µg/ml, Human
Immunohistochemistry(Paraffin-embedded Section), 2-5
µg/ml, Human Immunocytochemistry/Immunofluorescence,
5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10⁶ cells,
Human

Expiration Date


12 months from date of receipt.



Western blot
analysis of
RPA32/RPA2
usin...



IHC analysis
of
RPA32/RPA2
using anti-
RP...



IHC analysis
of
RPA32/RPA2
using anti-
RP...