

Product Datasheet

POR Antibody (monoclonal, 7F5) (orb763206)

Description

POR Antibody (monoclonal, 7F5)

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FC, IHC, WB

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human POR (633-668aa RNMARDVQNTFYDIVAELGAMEHAQAVDYIKKLMTK), different from the related mouse and rat sequences by five amino acids.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Mouse IgG2b

Clonality

Monoclonal

MW

77 kDa

Uniprot ID


P16435

Dilution Range


Western blot, 0.25-0.5 µg/ml, Human
Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Human

Expiration Date

12 months from date of receipt.



Western blot analysis of POR using anti-...



IHC analysis of POR using anti-POR antib...



IHC analysis of POR using anti-POR antib...