

Product Datasheet

Transcription factor E2F2 E2F2 Antibody (orb76323)

Description

Transcription factor E2F2 E2F2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human E2F2 (422-427aa ISDLFDSYDLGDLIN), identical to the related mouse sequence.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

47506 MW

Uniprot ID

[Q14209](#)


Dilution Range

Western blot, 0.1-0.5µg/ml, Human, Rat, Mouse

Expiration Date

12 months from date of receipt.

1 2 3 4 5 6 7 8 9



Western blot analysis of Rat Lung Tissue...