


Product Datasheet

Blood Group A1B antibody (orb763358)

Description	Mouse monoclonal antibody to Blood Group A1B
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	AGG, FC
Immunogen	Glycoprotein fraction isolated from saliva of secretors with blood group A.
Target	Blood Group A1B
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgM
Clonality	Monoclonal
Purity	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Dilution Range	FC: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Flow
cytometry
surface
staining
pattern ...