



Product Datasheet

mCherry antibody (orb763362)


Description	Rabbit polyclonal antibody to mCherry
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	mCherry protein from Anaplasma marginale
Target	mCherry
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml, extracellular staining or intracellular staining - depending on expression. Western blotting: Recommended dilution: 1-2 µg/ml.
Isotype	Rabbit polyclonal
Clonality	Monoclonal
Purity	Purified by ligand affinity chromatography.
Uniprot ID	X5DSL3
Dilution Range	FC: 1-4 µg/ml, WB: 1-2 µg/ml
Expiration Date	12 months from date of receipt.



Flow
cytometry
surface
staining
pattern ...



Separation of
HEK293T/17
cells co-
transf...



Western
blotting
analysis of
mCherry flu...