



Product Datasheet

C-Myc-Tag antibody (orb763372)


Description	Mouse monoclonal antibody to C-Myc-Tag
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB 1:5,000 IP 1:200 IF 1:1,000
Expiration Date	12 months from date of receipt.



IP antibody
use: 3ug Flag
Mouse IgG1
per...



1ug C-Myc
fusion protein
+ Primary
antib...



IF analysis of
293 cells
transfected
wit...