


Product Datasheet

mCherry-Tag antibody (orb763378)

Description	Mouse monoclonal antibody to mCherry-Tag
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB 1:5,000 IF 1:200
Expiration Date	12 months from date of receipt.



1ug mCherry
fusion
protein +
Primary ant...