


Product Datasheet

Strep-Tag antibody (orb763384)

Description	Mouse monoclonal antibody to Strep-Tag
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB 1:5,000
Expiration Date	12 months from date of receipt.



1ug Strep II
fusion
protein +
Primary an...