



Product Datasheet

alpha-tubulin antibody (orb763399)


Description	Mouse monoclonal antibody to alpha-tubulin
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	52kDa
Uniprot ID	P68366
Entrez	7277
Dilution Range	WB 1:5,000-10,000 IP 1:200 IF 1:100-200 IHC 1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) HeLa, 2) Rat...



Immunohistochemical analysis of paraffin...



1: Input: Mouse Brain Tissue Lysate 2: I...