



Product Datasheet

HSP70 antibody (orb763407)


Description	Mouse monoclonal antibody to HSP70
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	70kDa
Uniprot ID	P0DMV8
Entrez	3303
Dilution Range	WB 1:1,000-2,000 IF 1:100-200 IHC 1:200-500
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Western blot analysis of extracts from H...



IF analysis of HeLa with orb763407 (Left...