



Product Datasheet

HER-2 antibody (orb763410)

Description	Mouse monoclonal antibody to HER-2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	138kDa
Uniprot ID	P04626
Entrez	2064
Dilution Range	WB 1:2,000-4,000 IHC 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) Hela, 2) Mou...



IHC staining of human breast cancer tiss...