


Product Datasheet

AMACR antibody (orb763412)

Description	Mouse monoclonal antibody to AMACR
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	42kDa
Uniprot ID	Q9UHK6
Entrez	23600
Dilution Range	WB 1:1,000 IHC 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) HepG2, 2) Mo...



IHC staining of mouse prostate adenocarc...