



Product Datasheet

GFAP antibody (orb763423)


Description	Mouse monoclonal antibody to GFAP
Species/Host	Mouse
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	45kDa
Uniprot ID	P14136
Entrez	2670
Dilution Range	WB 1:2,000-5,000 IHC 1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat Brain Tissue...



Western blot analysis of Mouse Brain Tissue...



Immunohistochemical analysis of paraffin...