



Product Datasheet

CK17 antibody (orb763445)

Description	Mouse monoclonal antibody to CK17
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	48kDa
Uniprot ID	Q04695
Entrez	3872
Dilution Range	WB 1:1,000 IHC 1:200 IP 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) HeLa, 2) MCF...



1: Input: HeLa Cell Lysate 2: IP product...