


Product Datasheet

CD25 antibody (orb763456)

Description	Mouse monoclonal antibody to CD25
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P01589
Entrez	3559
Dilution Range	IHC 1:200-500
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...