



Product Datasheet

Cleaved PARP antibody (orb763472)

Description	Mouse monoclonal antibody to Cleaved PARP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	116,89kDa
Uniprot ID	P09874
Entrez	142
Dilution Range	WB 1:1,000-3,000 IHC 1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) Jurkat , 2) ...



Immunohistochemical analysis of paraffin...