



Product Datasheet

alpha skeletal muscle actin antibody (orb763481)


Description	Mouse monoclonal antibody to alpha skeletal muscle actin
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	42kDa
Uniprot ID	P68133
Entrez	58
Dilution Range	WB 1:5,000-10,000 IP 1:200 IHC 1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) Hela, 2) Mou...



1: Input: Mouse Brain Tissue Lysate 2: I...



Immunohistochemical analysis of paraffin...