



## Product Datasheet

### beta II tubulin antibody (orb763482)

|                            |                                                                                                            |
|----------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to beta II tubulin                                                               |
| <b>Species/Host</b>        | Mouse                                                                                                      |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                                               |
| <b>Tested Applications</b> | IP, WB                                                                                                     |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                          |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4.<br>Store at -20°C. Avoid repeated freeze-thaw cycles. |
| <b>Concentration</b>       | 1mg/ml                                                                                                     |
| <b>Storage</b>             | Store at -20°C. Avoid repeated freeze-thaw cycles.                                                         |
| <b>Note</b>                | For research use only                                                                                      |
| <b>Isotype</b>             | IgGq                                                                                                       |
| <b>Clonality</b>           | Monoclonal                                                                                                 |
| <b>MW</b>                  | 50kDa                                                                                                      |
| <b>Uniprot ID</b>          | <b>P68371</b>                                                                                              |
| <b>Entrez</b>              | <b>7280</b>                                                                                                |
| <b>Dilution Range</b>      | WB 1:50,000-100,000 IP 1:200                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                            |



Western blot analysis of 1) Mouse Brain



1: Input: Mouse Brain Tissue Lysate 2: I...