


Product Datasheet

Flotillin-1 antibody (orb763487)

Description	Mouse monoclonal antibody to Flotillin-1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	49kDa
Uniprot ID	075955
Entrez	10211
Dilution Range	WB 1:1,000-2,000
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) Mouse Brain

...