



Product Datasheet

GAP43 antibody (orb763510)

Description	Mouse monoclonal antibody to GAP43
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	38,43kDa
Uniprot ID	P17677
Entrez	2596
Dilution Range	WB 1:1,000-2,000 IHC 1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) Hela, 2) 293...



Immunohistochemical analysis of paraffin...