



Product Datasheet

GRP78 antibody (orb763533)


Description	Mouse monoclonal antibody to GRP78
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	78kDa
Uniprot ID	P11021
Entrez	3309
Dilution Range	WB 1:1,000-2,000 IHC 1:100-200
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) Hela, 2) Rat...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...