



Product Datasheet

Phospho-ERK 1/2 antibody (orb763591)

Description	Mouse monoclonal antibody to Phospho-ERK 1/2
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	44,42kDa
Uniprot ID	P27361 , P28482
Entrez	5595 , 5594
Dilution Range	IHC 1:100-200
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...