


Product Datasheet

GFP-Tag antibody (orb763629)

Description	Mouse polyclonal antibody to GFP-Tag
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB 1:5,000 IHC 1:200 IF 1:200 IP 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of GFP Recombinant...