


## Product Datasheet

### GAPDH antibody (orb763633)

|                            |                                                           |
|----------------------------|-----------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to GAPDH                       |
| <b>Species/Host</b>        | Rabbit                                                    |
| <b>Reactivity</b>          | Gallus, Human, Monkey, Mouse, Rabbit, Rat, Sheep, Xenopus |
| <b>Conjugation</b>         | Unconjugated                                              |
| <b>Tested Applications</b> | IF, IHC, WB                                               |
| <b>Immunogen</b>           | Recombinant Protein                                       |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4.      |
| <b>Concentration</b>       | 1mg/ml                                                    |
| <b>Storage</b>             | Store at -20°C. Avoid repeated freeze-thaw cycles.        |
| <b>Note</b>                | For research use only                                     |
| <b>Isotype</b>             | IgG                                                       |
| <b>Clonality</b>           | Polyclonal                                                |
| <b>MW</b>                  | 37kDa                                                     |
| <b>Uniprot ID</b>          | <b>P04406</b>                                             |
| <b>Dilution Range</b>      | WB 1:5,000 IHC 1:200                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                           |



Western blot analysis of 293T (1), Rat b...