


Product Datasheet

EFHD1 antibody (orb763639)

Description	Rabbit polyclonal antibody to EFHD1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	27kDa
Uniprot ID	Q9BUP0
Dilution Range	WB 1:1,000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T cell lysat...