


Product Datasheet

GAPDH antibody (orb763645)

Description	Rabbit polyclonal antibody to GAPDH
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37kDa
Uniprot ID	P04406
Dilution Range	WB 1:5,000-10,000 IHC 1:200 IF 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of Zebrafish skeletal tissue.