



Product Datasheet

Angiotensin (1-7) Mas Receptor antibody (orb763707)

Description	Rabbit polyclonal antibody to Angiotensin (1-7) Mas Receptor
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36-50kDa
Uniprot ID	P04201
Dilution Range	WB 1:2,000-5,000 IHC 1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis
of 1) Human Brain ...



Immunohistochemical
analysis of paraffin...