


Product Datasheet

TLR8 antibody (orb763715)

Description	Rabbit polyclonal antibody to TLR8
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NR97
Dilution Range	IHC 1:200-500
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...