

Product Datasheet

Human IgM (orb76373)

Description	Mouse monoclonal antibody to Human IgM
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	ChromPure Human IgM (myeloma), whole molecule.
Preservatives	Purified antibody in PBS containing 0.03% sodium azide.
Concentration	1mg/ml
Storage	4°C; -20°C for long term storage. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.