



## Product Datasheet

### GAPDH antibody (orb763749)

|                            |                                                      |
|----------------------------|------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to GAPDH                  |
| <b>Species/Host</b>        | Rabbit                                               |
| <b>Reactivity</b>          | Bacteria, Human, Mouse, Rat                          |
| <b>Conjugation</b>         | Unconjugated                                         |
| <b>Tested Applications</b> | IHC, WB                                              |
| <b>Immunogen</b>           | Recombinant Protein                                  |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. |
| <b>Concentration</b>       | 1mg/ml                                               |
| <b>Storage</b>             | Store at -20°C. Avoid repeated freeze-thaw cycles.   |
| <b>Note</b>                | For research use only                                |
| <b>Isotype</b>             | IgG                                                  |
| <b>Clonality</b>           | Polyclonal                                           |
| <b>MW</b>                  | 37kDa                                                |
| <b>Uniprot ID</b>          | <b>P04406</b>                                        |
| <b>Dilution Range</b>      | WB 1:5,000-10,000 IHC 1:200-500                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                      |



Western blot analysis of HeLa with GAPDH...



Immunohistochemical analysis of paraffin...