



Product Datasheet

STAT5a antibody (orb763752)

Description	Rabbit polyclonal antibody to STAT5a
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90kDa
Uniprot ID	P42229
Dilution Range	WB 1:500-1,000 IHC 1:100-200
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 with STAT5...



Immunohistochemical analysis of paraffin...