



Product Datasheet

Flotillin-2 antibody (orb763772)

Description	Rabbit polyclonal antibody to Flotillin-2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49,58kDa
Uniprot ID	Q14254
Dilution Range	WB 1:1,000-2,000 IHC 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) A431, 2) Mou...



Immunohistochemical analysis of paraffin...