


Product Datasheet

MEK4 (F194) antibody (orb763869)

Description	Rabbit polyclonal antibody to MEK4 (F194)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44kDa
Uniprot ID	P45985
Dilution Range	WB 1:1,000-2,000
Expiration Date	12 months from date of receipt.



Western blot analysis of
1) A431 Cell
Ly...