

Product Datasheet

Cleaved-Caspase-3 p12 (D175) antibody (orb763876)

Description	Rabbit polyclonal antibody to Cleaved-Caspase-3 p12 (D175)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human Caspase 3. AA range:157-206
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20+35
Expiration Date	12 months from date of receipt.