



Product Datasheet

Cleaved-C1r LC (I464) antibody (orb763888)

Description	Rabbit polyclonal antibody to Cleaved-C1r LC (I464)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human C1R. AA range:445-494
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	27
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
various cells
U...



Western
blot analysis
of lysates
from 29...