


## Product Datasheet

### Cleaved-Cathepsin A 32k (R326) antibody (orb763899)

|                            |                                                                                                  |
|----------------------------|--------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Cleaved-Cathepsin A 32k (R326)                                     |
| <b>Species/Host</b>        | Rabbit                                                                                           |
| <b>Reactivity</b>          | Human, Mouse                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                                     |
| <b>Tested Applications</b> | ELISA, WB                                                                                        |
| <b>Immunogen</b>           | The antiserum was produced against synthesized peptide derived from human PPGB. AA range:277-326 |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4.                                             |
| <b>Concentration</b>       | 1 mg/ml                                                                                          |
| <b>Storage</b>             | Store at -20°C. Avoid repeated freeze-thaw cycles.                                               |
| <b>Note</b>                | For research use only                                                                            |
| <b>Isotype</b>             | IgG                                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                                       |
| <b>MW</b>                  | 33                                                                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                  |



Western  
Blot  
analysis of  
various  
cells u...



Western  
blot  
analysis of  
lysates  
from 29...