


Product Datasheet

Cleaved-MPO 89k (A49) antibody (orb763945)

Description	Rabbit polyclonal antibody to Cleaved-MPO 89k (A49)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from Cleaved-MPO 89k (A49) . at AA range: 40-120
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	89
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
various cells
u...