



Product Datasheet

Ref-1 (Acetyl Lys6) antibody (orb763972)


Description	Rabbit polyclonal antibody to Ref-1 (Acetyl Lys6)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized Acetyl-peptide derived from human APE1 around the Acetylation site of Lys6. AA range:1-50
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35
Expiration Date	12 months from date of receipt.



Western Blot analysis of A549 cells using...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...