

Product Datasheet

Tubulin alpha (Acetyl Lys40) antibody (orb763989)

Description Rabbit polyclonal antibody to Tubulin alpha (Acetyl Lys40)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC-P, WB

Immunogen The antiserum was produced against synthesized Acetyl-peptide derived from human TUBA1A around the Acetylation site of Lys40. AA range:1-50

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 50


Expiration Date 12 months from date of receipt.



Western Blot analysis of Jurkat cells us...



Western Blot analysis of JK cells using ...



Immunohistochemical analysis of paraffin...