

Product Datasheet

Tubulin alpha (Acetyl Lys352) antibody (orb763990)

Description

Rabbit polyclonal antibody to Tubulin alpha (Acetyl

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, WB

Applications

Immunogen

The antiserum was produced against synthesized Acetyl-peptide derived from human TUBA1B around the Acetylation site of Lys352. AA range:311-360

Preservatives

PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration

1 mg/ml

Storage

Store at -20°C. Avoid repeated freeze-thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal


MW

50


Expiration Date

12 months from date of receipt.


KDa K562 / A549 / Jurkat



Western Blot analysis of K562, A549, Jur...



Western Blot analysis of K562 cells usin...



Western blot analysis of lysate from K56...