


Product Datasheet

Cytokeratin-pan (Acetyl Lys194) antibody (orb763992)

Description	Rabbit polyclonal antibody to Cytokeratin-pan (Acetyl Lys194)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized Acetyl-peptide derived from human Keratin-pan around the Acetylation site of Lys194. AA range:151-200
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	65
Expiration Date	12 months from date of receipt.

HeLa
117-
85-
75-
55-
45-
35-
25-
15-
5-
AcetylCytokeratin-pan (K194)

Western Blot analysis of HeLa cells using...



Western blot analysis of lysate from HeL...