


Product Datasheet

Tubulin alpha (Acetyl Lys112) antibody (orb763995)

Description	Rabbit polyclonal antibody to Tubulin alpha (Acetyl Lys112)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized acetyl-peptide derived from human Tubulin α around the acetylation site of K112.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50
Expiration Date	12 months from date of receipt.



Western
Blot
analysis of
HepG2 cells
usi...