



Product Datasheet

Ub (Acetyl Lys27) antibody (orb763998)


Description	Rabbit polyclonal antibody to Ub (Acetyl Lys27)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized acetyl-peptide derived from the human Ub around the acetylation site of K27.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Western Blot analysis of HeLa cells using Ub (Acetyl Lys27) antibody.



Western Blot analysis of HEPG2 cells using Ub (Acetyl Lys27) antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using Ub (Acetyl Lys27) antibody.