



Product Datasheet

NF-E4 (Acetyl Lys43) antibody (orb763999)


Description	Rabbit polyclonal antibody to NF-E4 (Acetyl Lys43)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized acetyl-peptide derived from the human NF-E4 around the acetylation site of K43.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	19
Expiration Date	12 months from date of receipt.



Western Blot analysis of AD-293 cells us...



Western Blot analysis of AD-293 cells us...



Immunohistochemical analysis of paraffin...