



Product Datasheet

PCAF (Acetyl Lys428) antibody (orb764002)


Description	Rabbit polyclonal antibody to PCAF (Acetyl Lys428)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized acetyl-peptide derived from the human PCAF around the acetylation site of K428.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	93
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 cells using PCAF (Acetyl Lys428) antibody.



Western Blot analysis of HepG2 cells using PCAF (Acetyl Lys428) antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using PCAF (Acetyl Lys428) antibody.