

Product Datasheet

Histone H1 antibody (orb764049)

Description	Rabbit polyclonal antibody to Histone H1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide of Histone H1
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17-25
Expiration Date	12 months from date of receipt.